

Laser Capture Microdissection

Pictures: Before, After, Cap

- . Dis-enable laser before taking pictures.
- . Re-focus image before each shot.

Spot Size: 7.5 μ m, 15 μ m, 30 μ m

- . Focus spot using lowest power magnification (10 X) first.
- . 7.5 μ m: .most effective for removal of smaller sections within the sample.
 - .Power: 50-55 mW
 - .Duration: 3.50- 4.00 ms
- . 15 μ m: .most effective for removal of epithelium.
 - .Power: 80-85 mW
 - .Duration: 1.5-2 ms
- .30 μ m: .Power: 45-55 mW
 - .Duration: 5.00 ms +

Sections: Frozen

- . Thin sections work best (~ 5 μ m thick, max 10 μ m).
- . Place section in center of slide.

Staining: LCM Kit, H&E, Quick Stain.

- . For best results, perform LCM immediately following staining.
- . Slides must be completely dry for effective LCM.
- . H&E staining yields least favorable results, while LCM kit staining (modified protocol) allows for optimal removal.
- . Re-Hydration of previously stained slides are not as effective for LCM.

LCM Staining: Protocol I

*Take precautions to ensure RNase free conditions

*Ensure that solutions are fresh (changed daily)

LCM Staining Kit (modified)

- 1) Remove sections from freezer, and allow them to thaw for 1 minute. If performing several rounds of staining at once, slides may be allowed to thaw in fumehood for up to 5 minutes.
 - 2) Place slides in 75% ethanol for 30 seconds.
 - 3) Wash with distilled water for 30 seconds.
 - 4) Using a staining rack, dip slides in Haematoxyllin solution for 15-20 seconds (~ 10 dips)
 - 5) Wash with distilled water for 30 seconds.
 - 6) Place slides in 75% ethanol for 30 seconds.
 - 7) Place slides in 95% ethanol for 30 seconds.
 - 8) Place slides in 100% ethanol for 30 seconds.
 - 9) In the fumehood, place slides in Xylene for 3 minutes.
 - 10) Dry slides in fumehood for 10 minutes. Gently remove any surrounding moisture with a kimwipe.
 - 11) Store slides in a slide box containing dessicant.
- ** Ensure that slides are completely dry before starting LCM.**

LCM Staining: Protocol II

***Take precautions to ensure RNase free conditions**

H & E staining

- 1) Place slides in 100% ethanol for 4 seconds.
- 2) Transfer slides to 90% ethanol solution for 4 seconds.
- 3) Transfer slides to 70% ethanol for 4 seconds.
- 4) Wash slides with distilled water
- 5) Place slides in Haematoxylin for 60 seconds.
- 6) Wash slides with distilled water.
- 7) Place slides in Lithium Carbonate for 10 seconds
- 8) Wash slides with distilled water
- 9) Place slides in Eosin for 20 seconds.
- 10) Place slides in 70% ethanol for 10 seconds.
- 11) Place slides in 90% ethanol for 10 seconds.
- 12) Place slides in 100% ethanol for 10 seconds.
- 13) Place slides in 100% ethanol for 10 seconds.
- 14) In fumehood, directly transfer slides to Xylene for 3 minutes.
- 15) Leave slides to dry in the fumehood for 10 minutes.
- 16) Place slides in a slide box with fresh dessicant.

**** Ensure that slides are completely dry before starting LCM.**

LCM Staining: Protocol III

*Take precautions to ensure RNase free conditions.

Quick Staining

- 1) Place slides in Fixative solution for 5 seconds (or 5 dips).
- 2) Place slides in Eosin for 5 seconds.
- 3) Place slides in Methylene blue for 5 seconds.
- 4) Wash slides under the water tap, on the opposite side of staining.
- 5) Allow slides to dry in fumehood for 10 minutes.
- 6) Place slides in slide box containing fresh dessicator.

**Ensure that slides are completely dry before starting LCM.

LCM Re-Hydration: Protocol

- 1) Place slides in 95% ethanol for for 60 seconds.**
- 2) Place slides in 100% ethanol for 30 seconds.**
- 3) In the fumehood place cells in Xylene for 3 minutes.**
- 4) Allow slides to dry for 10 minutes in fumehood.**

LCM- Pre-mRNA Isolation: Protocol

- 1) Place cap into 500 μ l microcentrifuge tube with 100 μ l Buffer RLT.
- 2) Invert tube several times to digest material off cap. (Alternative: incubate inverted at 42 C for 30 minutes).
- 3) Vortex cap end of tube.
- 4) Centrifuge at 2000 x g for 60 seconds.
- 5) Remove cap.
- 6) Adjust the volume to 350 μ l with Buffer RLT.
- 7) Add one volume of 70 % ethanol.
- 8) Proceed using Rneasy Mini kit for isolation of total mRNA from animal cells.

To further concentrate the Sample (optional)

- 1) Concentrate the sample with 0.5 ammonium acetate and 5 μ g glucogen.
- 2) Centrifuge
- 3) Wash pellet with 75% and 100% ethanol.
- 4) Allow to air dry, briefly.
- 5) Solubilize in 10 μ l Rnase-free water.